


Viewing Your Individual Dose Report

Issue Date: 02/22/18

Revision Date: 02/22/18

1. Go to Landauer’s website: <https://myLDR.com>.

Response: The following Web page is displayed:



2. Enter the login information as per your area of work and account:

- Username:
 - Clinical
 - UM207111 (Michigan Medicine - Monthly Accounts)
 - UM207112 (Michigan Medicine - Quarterly Accounts)
 - Research
 - UM102200 (Campus – Monthly Accounts)
 - UM102400 (Campus – Quarterly Accounts)
- Password: MyDose1

3. Enter the **Account Number** (without the “UM”), then enter the **Serial Number** on your dosimeter. The illustration below shows where the serial number is located on each type of dosimeter.



Response: Dose information will appear as follows without any identifying data. This summary provides your most recent quarterly results as well as the year to date total for the current year and the total for the previous year.

To protect your privacy, no personal information is displayed

History Results
The doses are displayed in mrem

Dose Results

Total Records: 4
Doses as of 2014/06/24 13:14 CST

Dose Period	Total DDE	Total LDE	Total SDE	Beta	Total Neutron	Extremity
Q22014	13	13	12			
2014	16	15	14			
2013	22	22	24			
Lifetime	361	416	518	51		M

Page 1 of 1 Prev Next

Export Result to: Excel Export View Details

- For more detailed information, click **View Details** to see a breakdown of your dose by wear period.

NOTE: The Lifetime Dose listed is only for dosimeters received under that specific account. Please contact Radiation Safety Services/EHS if you would like to obtain your complete Lifetime Dose History.

History Details

Dosimeter Dose

Total Records: 19
Doses as of 2014/06/24 13:16 CST

Account	Subaccount Name	Subaccount Code	Begin Wear Date	End Wear Date	Dosimeter Type	Dosimeter Location	Total DDE	Total LDE	Total SDE	Beta	Total Neutron	Fast Neutron	Thermal Neutron	Extremity	Frequency	Control Subtracted	Not in Assign Dose
88509	HPO	HPO	2014/07/01	2014/07/31	Pa	Chest									Monthly	<input type="checkbox"/>	<input checked="" type="checkbox"/>
88509	HPO	HPO	2014/06/01	2014/06/30	Pa	Chest									Monthly	<input type="checkbox"/>	<input checked="" type="checkbox"/>
88509	HPO	HPO	2014/05/01	2014/05/31	Pa	Chest	11	11	10						Monthly	<input checked="" type="checkbox"/>	<input type="checkbox"/>
88509	HPO	HPO	2014/04/01	2014/04/30	Pa	Chest	2	2	2						Monthly	<input checked="" type="checkbox"/>	<input type="checkbox"/>
88509	HPO	HPO	2014/03/01	2014/03/31	Pa	Chest	1	1	1						Monthly	<input checked="" type="checkbox"/>	<input type="checkbox"/>
88509	HPO	HPO	2014/02/01	2014/02/28	Pa	Chest	2	1	1						Monthly	<input checked="" type="checkbox"/>	<input type="checkbox"/>
88509	HPO	HPO	2014/01/01	2014/01/31	Pa	Chest	M	M	M						Monthly	<input checked="" type="checkbox"/>	<input type="checkbox"/>
88509	HPO	HPO	2013/12/01	2013/12/31	Pa	Chest	M	M	M						Monthly	<input checked="" type="checkbox"/>	<input type="checkbox"/>

Page 1 of 1 Prev Next

Export Result to: Excel Export

- At the top right corner of the Web page, click **End Session**.